

Drug Dipstick (urine) (Amphetamine / Barbiturates / Benzodiazepines / Cocaine /Ketamine / Morphine / Marijuana / Tramadol)


DIRECTIONS FOR USE

1- Bring the pouch to room temperature before opening it. Remove the Test Dipstick from the sealed pouch and use it within one hour.

2- With arrows pointing toward the urine specimen, **immerse the test dipstick vertically in the urine specimen for at least 10-15 seconds**

3- Place the Test Dipstick on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear.

Read results at 5 minutes. Do not interpret the result after 10 minutes.



INTERPRETATION OF RESULTS


NEGATIVE: * **Two lines appear.** One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T)

***NOTE:** The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: **One colored line appears in the control line region (C).** No line appears in the test line region (T). This positive result indicates that concentration exceeds the detectable level

INVALID: **Control line fails to appear.** Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure.

Review the procedure and repeat the test using a new test.



Cut off concentration

No.	Item	Cut off concentration
1	Amphetamine (AMP)	1000 ng/ml
2	Barbiturates (BAR)	300 ng/ml
3	Benzodiazepines (BZO)	300 ng/ml
4	Cocaine (COC)	300 ng/ml
5	Ketamine (KET)	1000 ng/ml
6	Morphine (MOP)	300 ng/ml
7	Marijuana (THC)	50 ng/ml
8	Tramadol (TML)	100 ng/ml

LIMITATION

1. Rapid Test Dipstick (Urine) provides only a qualitative, preliminary analytical result. A secondary (GC/MS) is the preferred confirmatory method.
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.

PRECAUTIONS

1. For medical and other professional in vitro diagnostic use only. Does not use after the expiration date.
2. The test should remain in the sealed pouch until use.
3. All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
4. The used test should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Does not use beyond the expiration date.



Diagnostics

Diagnostic kits & Biological products