


## COVID-19 IgG/IgM Rapid Test (Whole Blood/Serum/Plasma)


### Whole Blood specimen:

- Transfer 20  $\mu$ l of specimen to the specimen well (S) with a pipette or a dropper
- add 2 drops of buffer (approximately 80  $\mu$ l) to Buffer well (B)




### For Serum or Plasma specimen

- Transfer 10  $\mu$ l of specimen to the specimen well (S) with a pipette or a dropper
- add 2 drops of buffer (approximately 80  $\mu$ l) to Buffer well (B)



### Result Interpretation

Wait for colored line(s) to appear. Read the result at 15 minutes. Don't interpret the result after 20 minutes.



### Precautions

- This package insert must be read completely before performing the test. Failure to follow directions in insert may yield inaccurate test results.
- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- Wear protective clothing such as laboratory coats, disposable gloves, and eye protection when specimens are assayed.
- Please ensure that an appropriate quantity of sample volumes are used for testing. Too much or too little sample size may lead to deviation of results.
- The used test should be discarded according to local regulations.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use.

**DO NOT FREEZE.** Do not use beyond the expiration date.